

### HEPG2 Subculturing in 100mm dish

1. Remove and discard culture medium.
2. Rinse cells with 3.5ml PBS and remove after 30 seconds to 1 minute
3. Add 1.5ml of Trypsin-EDTA solution to dish and place cells at 37 degrees to facilitate dissociation. Observe cells under an inverted microscope until cell layer is dispersed (usually within 3 to 8 minutes).
4. Add 3.5 mL of complete growth medium and collect cells in a 50ml conical tube.
5. Aseptically syringe cells 3x using a 5ml syringe and 22 gauge needle. **Discard both needle and syringe in sharps container. Do not disconnect needle from syringe, Do not recap needle.**
6. Add fresh complete growth media to new unused 100mm dishes
7. Add appropriate aliquots of the cell suspension to new dishes. (Recommended to use 1/2, 1/3, 1/4 split ratios)
8. Place cells in 37 degree incubator and shake in cross pattern to evenly distribute cells.

### Preparing gelatin coated dishes for HEPG2

1. Retrieve and open new unused culture dishes
2. Add only enough 0.1% gelatin to cover the surface of the dish (100mm use ~3mls, 6 wells use ~0.5-1ml/well)
3. Place the gelatin coated dishes in a 37 degree incubator for 10 minutes
4. Aspirate the gelatin.
5. The dishes are now coated and ready for media/cells.